## **Ultra Centrifuge**

**Make and Model**: Beckman Coulter (Optima XPN100)

**Applications**: The ultracentrifuge employs high-speed rotation to achieve significant g-force separations. Operating within the range of 100,000 xg to 1,000,000 xg, ultracentrifuges are instrumental in purifying and characterizing a spectrum of substances, ranging from low-molecular-weight polymers to multi-mega Dalton protein complexes and organelles.

**Charges**: For more information and special offers for students, MSME, SSI, National R&D Labs, please visit our website link:

https://www.niist.res.in/testing-and-analysis